





	No	Type		No	Type		No	Type
<b>Dental Ratio</b>								
RP, WP Tissue Level	55014	4						
NC Bone Level	55061	53						
RC Bone Level	55062	54						
			<b>Futura implants</b>			<b>IML Swiss Dental</b>		
			PRE 4.0	55014	4	Starfly 3.5-6.0	55034	26
			PRI 4.0	55037	29	Universe 2.9	55039	31
						Universe 3.4-5.0	55024	15
<b>Dental Tech</b>								
FT3 Implastic 3.25	55039	31						
FT3 Implastic 3.75	55049	41	<b>galimplant</b>					
FTP Implastic	55019	9	ICI 3.2 & 3.8	55046	38	<b>Implance</b>		
HX3 Implastic	55011	1	iPx 4012	55046	38	Tissue level	55014	4
TR2 Implastic	55014	4	MUSA	55128	94	Regular Bone level	55037	29
TRW2 Implastic	55014	4				Bone level 3.3	55065	57
			<b>GC Implant</b>					
<b>Dentalis</b>			Narrow Aadva	55081	73	<b>Implant Direct</b>		
Regular	55051	43	Regular Aadva	55066	58	Legacy 1/2/3/4, 3.0/3.2	55012	2
			Wide Aadva	55044	36	Legacy 1/2/3/4, 3.5/3.7-5.7/7.0	55051	43
			Short Aadva 4.2	55044	36	RePlant 3.5	55022	12
<b>Dentegris</b>			Plus Genesio 3.4-5.0	55085	77	RePlant 4.3-5.0	55010	13
3.3	55034	26	Plus SETIO 3.8	55011	1	RePlant 6.0	55032	24
3.75	55057	49				ReActive 3.5	55022	12
4.5	55049	41	<b>Geass</b>			ReActive 4.3-5.0	55010	13
3.8/4.1	55049	41	Way Extra	55081	73	3.0 Interactive	55046	38
5.5	55020	10				3.4 Interactive	55037	29
			<b>Global D</b>			3.0 Simply Iconic	55046	38
Cleanlant S-Clean 3.7-4.8	55020	10	In-Kone 4.0-5.0	55024	15	3.4 Simply Iconic	55037	29
Narrow S-Clean OneQ 3.3	55064	56	TwinKon 3.5-4.5	55122	88	Screwindirect 5.0 PF 3.7-5.7	55070	62
<b>Dentium</b>			<b>gmi</b>			<b>Implant One</b>		
Superline	55037	29	RP avantgard 3.75-4.25	55066	58	Std Thread 100 Series	55029	21
NR Line	55066	58	WP avantgard 4.75-5.75	55075	67	Std Thread 200 Series	55029	21
Simple Line II GH 1, 6	55014	4	RP, WP frontier	55051	43	Std Thread (blue) 300 Series	55029	21
			NP phoenix 3.3	55038	30	Std Thread (purple) 400 Series	55035	27
			Std-wide phoenix 3.3-5.0	55011	1	Std Thread (green) 500 Series	55035	27
<b>Dentsply</b>			monolith 3.0	55043	35	Wide Thread (blue) 300 Series	55029	21
PrimeTaper/OmniTaper 3.0	55058	50				Wide Thread (purple) 400 Series	55035	27
PrimeTaper/OmniTaper 3.6	55046	38	<b>GP Implants</b>			Wide Thread (green) 500 Series non-EU markets only	55129	95
PrimeTaper/OmniTaper 4.2	55057	49	SSI Spiral Shape Impant	55051	43			
PrimeTaper/OmniTaper 4.8	55055	47				<b>Implantswiss</b>		
PrimeTaper/OmniTaper 5.4	55055	47				Tissue level	55014	4
			<b>gt medical</b>			Regular Bone Level	55037	29
<b>DEVA</b>			HXE	55011	1	Bone Level 3.3	55065	57
Bone Level 3.5	55064	56	RP HXE	55012	2	Bone Level 2.9	55058	50
Bone Level 4.5	55051	43	WP HXE	55013	3			
			HXI 4.1	55037	29			
			OO2 HXI 4.7	55049	41	<b>Implantvision</b>		
<b>DIO</b>			OCT 4.3	55014	4	Internal Hex 3.7-5.7	55034	26
Submerged SM 3.8	55049	41						
Submerged SM 5.0	55020	10	<b>Hahn</b>					
Narrow UF II	55064	56	Tapered 3.0	55058	50	<b>Innobiosurg</b>		
Regular UF II	55020	10	Tapered 3.5	55046	39	Regular MagiCore	55018	8
			Tapered 4.3	55046	38			
<b>Ditron Dental</b>			Tapered 5.0	55037	29	<b>Intoss Anchor</b>		
MPI 3.3-6.0	55087	79	Tapered 7.0	55048	40	Internal hex 3.7-5.7	55088	80
Ultimate 3.75-6.0	55087	79						
						<b>Intra-lock</b>		
<b>DSI</b>			<b>Hi-Tec</b>			Intrahex 3.4 - 4.5	55051	43
Regular Internal Hex	55049	41	NTI, SLT Tapered Self Thread 3.3, 3.75	55034	26	Gold & Blue 3.4 - 6.0	55051	43
Narrow Internal Hex	55039	31	GTI, WTI Tapered Self Thread	55051	43	IndEx 3.4 - 3.75	55015	5
Regular Conical connection	55037	29	TPI Spark	55034	26	IndEx 4.0 - 4.3	55020	10
Narrow Conical connection	55046	38	Plus Logic	55034	26	FlatTop 3.75	55012	2
			LGI Logic 3.0-3.5	55046	38			
			LGI Logic 4.3-5.0	55037	29	<b>IRES</b>		
<b>DSP Biomedical</b>						Narrow External Hex 3.3	55012	2
External Standard (HE) Hexagonal	55014	4	<b>Hiossen</b>			Regular External Hex 3.7-5.2	55011	1
Internal Standard (HIS) Hexagonal	55034	26	HS system	55028	20	Narrow Internal Hex 3.3	55087	79
Internal Morse (HIM) Hexagonal	55046	38	mini ET system 3.2	55047	39	Regular Internal Hex 3.7-5.2	55049	41
Internal Large (HIL) Hexagonal	55037	29	mini ET system 3.5	55015	5	Regular Internal Octagon 3.7-5.2	55014	4
Morse Taper (CMI) Indexed	55025	16	regular ET system 4.0-7.0	55037	29			
			single EK system 3.5	55046	38	<b>Izen Implants</b>		
<b>DYNA Dental</b>			single EK system 4.0-7.0	55044	36	Narrow I-System X	55064	56
Helix DC	26	55034				Narrow I-System	55047	39
Helix ST	26	55034	<b>IDI Evolution</b>			Mini/Regular/Wide I-System	55020	10
Helix TM	16	55025	Standard Implant	55004	18	Narrow T-System	55047	39
			High B1One	55080	72	Mini T-System	55015	5
<b>Easy Implant</b>			Low B1One	55033	25	Regular/Wide T-System	55037	29
Ocean	55046	38				Mini/Regular/Wide R-System	55051	43
Lilas	55037	29	<b>IDI Implants</b>					
			BIO 3.5-5.0	55037	29	<b>J Dental Care</b>		
<b>Edierre</b>			ST ID Cam 4.2	55010	13	JD Evolution 3.7-6.0	55067	59
3-lobe	55021	11	All ID	55010	13	Plus+ JD Evolution 3.7-6.0	55067	59
4-lobe	55025	16	ID3 3.1	55021	11	S JD Evolution 3.25	55044	36
						S JD ICON ultra 2.75	55064	56
<b>Edison Medical</b>			<b>ids/Biotech dental</b>			Conical Abutment	55080	72
MUA	55033	25	Contact, Contact-S 3.6 - 4.8	55045	37			
			Aspire Reflect	55081	73	<b>Kamal Medtech</b>		
<b>Euroteknika</b>			Certus Reflect	55019	9	i-fix 3,5-5,0	55037	29
Naturactis 3.5-5.0	55081	73	Tapered Reflect	55034	26	s-fix 3,6-6,0	55046	38
Naturactis/Naturall+ 3.0	55064	56	Recover Reflect	55020	10			
Uneva	55011	1	Rapid Reflect	55020	10			



	No	Type		No	Type		No	Typ
<b>Keystone</b>								
I-HEXMRT 7.0-9.0	55051	43						
Narrow Prima/Connex	55049	41						
Regular Prima/Connex	55067	59	<b>Meisinger Implants</b>			<b>Novem</b>		
Wide Prima/Connex	55041	33	RP, WP OKTAGON Tissue level	55014	4	3.4	55039	31
Narrow Prima/Genesis	55049	41	NC OKTAGON Bone Level	55061	53	3.8	55016	6
Regular Prima/Genesis	55067	59	RC OKTAGON Bone Level	55062	54	4.2, 5.0	55057	49
Wide Prima/Genesis	55041	33						
Narrow Restore	55079	71	<b>Microdent</b>			<b>NucleOss</b>		
Regular Restore	55011	1	Universal	55011	1	Tpure Bone Level	55004	18
Wide Retore	55013	3	Microdent 2.8-3.25	55068	60	T6 Bone Level	55074	66
TILOBEMAXX 7.0-9.0	55041	33	Microdent 3.3-5.0	55011	1			
			Microdent 5.5-8.0	55013	3			
			Genius 3.0	55025	16	<b>OCO Biomedical</b>		
<b>Klockner</b>			Genius 3.5-5.0	55051	43	Regular Engage 3.25/3.0	55027	19
Micro VEGA 3.0	55074	66	Ektos 3.2	55023	14	Regular Engage 4.0/3.5	55051	43
Narrow VEGA 3.5	55065	57	Ektos 3.7-8.0	55051	43	Regular Engage 5.0/4.5	55051	43
Regular VEGA 4.0, 4.5	55034	26	Trylogic 3.5-4.2	55051	43	Regular Engage 6.0/5.7	55051	43
Narrow KL 3.3	55012	2	Trylogic 5.0	55075	67	Engage Plus 4.0-6.0/3.25		
Regular KL 3.7, 4.2	55011	1				ERI 3.25/3.7	55069	61
Wide KL 4.7	55013	3	<b>MIS</b>			ERI 4.0/4.7	55014	4
SK2	55038	30	Narrow (C1, V3) Conical Connection	55081	73	ERI 5.0/5.7	55014	4
Essential Cone	55067	59	(C1, V3) Conical Connection	55087	79	TSI 3.25/3.7	55068	60
S3M 3.1	55079	71	(C1) Conical Connection	55041	33	TSI 4.0/4.7	55011	1
			Narrow Seven and M4	55027	19	TSI 5.0/5.7	55011	1
<b>Konmet</b>			Standard Seven and M4	55049	41			
Classic 4.0, 4.8	55014	4	Wide Seven 5.0	55051	43	<b>OsseoFuse</b>		
Cylindrical 3.5-4.0	55065	57	NP, SP, WP Connect Abutment	55115	81	Hexapulus/ HexaPlus-S Sinus 3.75-7.5	55081	73
Cylindrical 4.5-5.5	55057	49	Multi Unit System	55033	25	ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5
Short Cylindrical 5.0-6.0	55078	70						
Conical 4.4/3.4	55065	57	<b>Mode Implant</b>			<b>Osstem</b>		
Conical 4.8/3.8, 5.5/4.5	55057	49	NP Rapid 3,3	55046	38	Regular US	55011	1
			NP Rapid 3,7	55046	38	Mini US	55012	2
<b>Kyocera</b>			RP Rapid 4,1	55027	19	Wide US	55013	3
NP FINESIA Bone Level 3.2-3.4	55058	50	RP Rapid 4,7	55027	19	Wide PS US	55011	1
RP FINESIA Bone Level 3.7-4.2	55034	26	RP Rapid 5,2	55027	19	Regular, Wide SS	55028	20
WP FINESIA Bone Level 4.7-5.2	55014	4	NP Shorter 3,7	55046	38	Regular, Regular Short TS 6, 7	55037	29
RP FINESIA Tissue Level 3.7-4.7	55014	4	RP Shorter 4,1	55027	19	Mini TS 3.5	55015	5
WP FINESIA Tissue Level 4.7-5.2	55014	4	RP Shorter 4,7	55027	19	Mini TS 3.0	55047	39
TP/ST POI (HAC) EX 3.7	55085	77	RP Shorter 5,2	55027	19			
TP/ST POI (HAC) EX 4.2	55025	16	NP Level 3,3	55027	19	<b>Osteo Plus</b>		
POI (HAC) 2P 3.7/4.2	55004	18	NP Level 3,7	55027	19	PF Iridium 3.8, 4.2	55079	71
			RP Level 4,1	55027	19	She 3.45, 3.75	55011	1
<b>Lasak</b>			RP Level 4,7	55027	19	Shi 3.75, 4.2	55049	41
QR BioniQ	55066	58	RP Level 5,2	55027	19	Shi 5.0	55051	43
QN BioniQ	55021	11	NP Short 3,7	55065	57			
			RP Short 4,1	55081	73	<b>OT Medical</b>		
<b>Leader</b>			RP Short 4,7	55081	73	OT-F2 3.4, 3.8	55053	45
Internal Hex 3.3-5.5	55034	26	RP Short 5,2	55081	73	OT-F2 4.1	55073	65
Cono-In	55067	59						
			<b>myplant GmbH</b>			<b>Overmed</b>		
<b>Maxillent</b>			myplant two 3.5-4.5	55077	69	Easyline Easydip	55020	10
iRaise	55049	41	myplant bio 4.0	55077	69	ETLINE	55034	26
<b>Medentika</b>			<b>Neo Biotech</b>			<b>Oxy Implant</b>		
NI Microcone	55021	11	IS II 4.0	55020	10	Internal	55020	10
RI Microcone	55081	73	IS III 4.0	55020	10	Narrow Kone Conical	55021	11
RI Quattrocone	55081	73				Regular Kone Conical	55051	43
30 AI Quattrocone	55081	73	<b>Neodent</b>			Regular PSK Line	55051	43
Procone 3.3	55119	85	Morse Connection Cone	55025	16	Mini External	55079	71
Procone 3.8	55031	23	WS Connection Titamax	55056	48	Regular External	55011	1
Procone 4.3	55031	23	Morse Connection Grand	55066	58	Maxi External	55011	1
Procone 5.0	55032	24	GM Micro Conical Abutment	55080	72			
			Helix Short Purple	55056	48	<b>P-I Brånemark</b>		
<b>Medentis</b>						HEX 3.3	55115	81
BL ICX Premium/Active Master 3.75-4.8	55062	54	<b>Neoss</b>			HEX 3.75-5.0	55117	83
Tissue Level ICX Plus 3.45	55039	31	SP 3.5 implants and wider	55004	18	AMP 3.3	55039	31
BL ICX Premium/Active Master 3.3	55058	50	NP 3.25	55027	19	AMP 3.75-4.8	55118	84
			Access Abutment	55060	52	MT 3.3-4.8	55046	38
						Conical Abutment HEX, AMP MT	55080	72
<b>Medical Instinct</b>			<b>Nobel Biocare</b>			<b>Paltop</b>		
Hex 3.0	55071	63	RP Brånemark System Ext Hex	55011	1	Standard 3.75-5.0	55034	26
Hex 3.4, 4.0, 5.0	55052	44	NP Brånemark System Ext Hex	55012	2	Narrow	55027	19
Cone 3.0, 3.4, 4.0	55071	63	WP Brånemark System Ext Hex	55013	3	Conical PCA	55066	58
Cone 5.0	55072	64	NP Tri-Lobe Connection	55022	12			
			RP, WP Tri-Lobe Connection	55010	13	<b>Park Dental Research</b>		
<b>Medigma Biomedical</b>			Tri-Lobe Connection 6.0	55032	24	P series Starvent	55046	38
Mars 3.5-6.0	55057	49	Conical Connection 3.0	55058	50	XP series Starvent	55049	41
FixTite 3.5-6.0	55057	49	Narrow Conical Connection	55046	38	Narrow ridge Startarius	55046	38
Fix-a-Dent 3.75-6.0			Regular Conical Connection	55037	29			
			Wide Conical Connection 5.5	55048	40	<b>Phibo</b>		
<b>Medimecca</b>			Regular N1	55126	92	NP Aurea EVO	55027	19
PT SLA Chaorum	55047	39	Narrow N1	55127	93	RP, WP Aurea EVO	55051	43
			NP, RP MUA	55080	72	TSA S3	55074	66
<b>Megagen</b>			NP, RP, WP On1	55042	34	TSA S4, S5	55034	26
AnyRidge	55051	43				TSH S3, S4, S5	55011	1
AnyOne Internal	55037	29	<b>Noris Medical</b>					
Mini	55064	56	Tuff 3.3-6.0	55034	26			
MUA	55080	72						
Octa	55121	87						



	No	Type		No	Type		No	Typ
<b>PLATON JAPAN</b>								
NC Eight-Lobe Pro 3.5-4.1	55046	38						
RC Eight-Lobe Pro 4.5-5.0	55014	4						
Type IV PRO	55020	10						
<b>Prodent Italia</b>			<b>Stark Implants</b>			<b>Xgate</b>		
Orange Twinner & Prime, Conomet, Collar	55027	19	Narrow Gravis 3.0-3.5	55064	56	Conical 3.3-5.0	55034	26
Fuchsia Twinner & Prime, Conomet, Collar	55051	43	RP Gravis 3.8-6.0	55037	29			
Green Twinner & Prime, Conomet, Collar	55067	59				<b>Xive</b>		
Yellow Twinner & Prime, Conomet, Collar	55020	10	<b>Straumann</b>			S-Plus 3.0	55030	22
Blue Prime, Conomet, Collar	55020	10	Regular & Wide neck Tissue level	55014	4	S-Plus 3.4	55061	53
			Narrow Neck Crossfit (NNC) Tissue level	55076	68	S-Plus 3.8-5.5	55062	54
			Tapered Narrow (NC) Bone level	55061	53	TG Plus 3.4-4.5	55018	8
			Tapered Regular (RC) Bone level	55062	54			
			Tapered SC Bone level 2.9	55064	56	<b>Zeramex</b>		
<b>Qualia</b>			BLX all platforms Bone level	55027	19	XT Implant 4.2, 5.5	55086	78
NeoActive	55046	38	TLX all platforms Tissue Level X	55115	81			
			NC, RC Screw retained abutment	55080	72	<b>Ziacom</b>		
<b>Quest</b>			PURE ceramic 4.1	55047	39	NP Internal hexagonal connection	55027	19
NP Hylock	55030	22				RP Internal hexagonal connection	55057	49
RP Hylock	55066	58	<b>Surgikor</b>			WP Internal hexagonal connection	55057	49
			Fixation 3.0	55058	50	NP External hexagonal connection	55079	71
<b>Resista</b>			Fixation 3.5	55066	58	RP External hexagonal connection	55117	83
WIS 3.0-5.0	55049	41	Fixation 4.3	55037	29	WP External hexagonal connection	55014	4
			Solution 3.5, 4.0	55066	58	RP Conical connection	55046	38
<b>Rex Implants</b>			Solution 4.5, 6.0	55037	29	WP Conical connection	55075	67
TL External hex 1.8	55060	52	Immediate all	55049	41			
TL Internal Hex 2.9	55034	26	Versatile all	55049	41	<b>Zimmer/Zimvie</b>		
BL Internal Hex 2,9	55034	26				Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
			<b>Sweden &amp;#03B; Martina</b>			Screw-Vent, Tapered Screw-Vent 4.7-6.0	55051	43
<b>Ritter</b>			Syra 3.8-5.0	55011	1	SwissPlus 3.8	55070	62
Multi Unit	55033	25	Outlink2 3.75-5.0	55011	1	SwissPlus 4.8	55014	4
Narrow Spiral SB/LA 2.9	55081	73	Outlink2 3.3	55079	71	Spline Twist 3.75-5.0	55038	30
Standard Spiral SB/LA 3.75	55049	41	Shelta Power 3.8-5.0	55025	16	NP Ezthetic 3.1	55065	57
			Shelta 3.8-5.0	55025	16	TSX 3,5	55049	41
<b>RS Implant</b>			Premium POWER 3.3-5.0	55025	16	TSX 4,5	55051	43
3.5-6.0	55020	10	Premium ONE 3.3-5.0	55025	16			
			Prama 3.3-5.0	55079	71	<b>Zuga</b>		
<b>SchütZ Dental</b>			CSR-N DAT	55023	14	All	55027	19
Impla cylindrical	55052	44	CSR, DAT connection	55035	27			
IMPLA micro Retention cone 4.2	55052	44	Shelta Abutment non-EU markets only	55128	94			
			Youse I	55034	26			
<b>SGS Implant</b>			Youse V-N	55065	57			
P1, P7	55049	41	Youse V-R	55020	10			
P1D, P5D, P7D	55049	41						
P7N	55059	51	<b>Tag Dental</b>					
			Axis	55051	43			
<b>SHOFU</b>			massif	55051	43			
BIOFIX 3.0	55064	56						
BIOFIX 4.0	55046	38	<b>Tatum</b>					
BIOFIX 4.5	55046	38	Integrity 3.7-5.0	55015	5			
			Integrity 6.0-8.0	55040	32			
<b>SIC Invent</b>			T Implant	55063	55			
MUA	55120	86						
SiCace 3.4-5.0	55052	44	<b>TBR</b>					
SiCmax, SiCtapered 3.7-5.2	55052	44	Infinity PF 3.5, 4.0	55004	18			
SiCvantage max, SiCvantage tapered 3.0	55026	17	M Implant	55044	36			
SiCvantage® max, SiCvantage tapered 3.7	55116	82						
SiCvantage max, SiCvantage tapered 4.2-5.2	55054	46	<b>Thommen</b>					
			SPI 3.5	55021	11			
<b>Simply Implants</b>			SPI 4.0	55065	57			
Internal Hex	55051	43	SPI 4.5-6.0	55066	58			
<b>SIN</b>			<b>TRATE</b>					
SW CM Strong	55023	14	All diam ROOTFORM	55065	57			
SW Int Hex Strong	55035	27						
Tryon 3.25-5.0	55011	1	<b>TRI Dental implants</b>					
CM Unittite 2.9	55048	76	Tissue level Octa	55014	4			
CM Unittite 3.5	55077	69	all Tri-Vent	55034	26			
CM Compact Unittite 4.0	55085	77	Narrow	55025	16			
			<b>Trinon</b>					
<b>Southern Implants</b>			Q2	55011	1			
Tri-Nex 4.3-5.0	55010	13	QK	55014	4			
Tri-Nex 3.5	55022	12	QZA 3.5	55049	41			
External Hex	55011	1	QZA 4.5, 5.7	55051	43			
MAX Ext Hex	55011	1						
0, 12, 24, 36 Co-Axis Ext Hex	55011	1	<b>URIS</b>					
Piccolo Ext Hex	55050	42	Narrow Omni 3.0-3.5	55115	81			
IT-Connection	55014	4	Regular Omni 4.0-6.5	55020	10			
Deep Conical 3.0	55064	56						
Deep Conical 3.5-4.0	55046	38	<b>Vulkan</b>					
Tapered Deep Conical 5.0	55017	7	Internal hex	55034	26			
			NP External hex	55012	2			
<b>Spiral Tech</b>			RP External hex	55011	1			
Hex connection ESI	55051	43	MIP Conical con	55058	50			
Hex connection Ultimate	55051	43	NP Conical con	55046	38			
Conical ESI 3.5	55044	36						
Conical ESI 4.3, 5.0	55037	29	<b>Winsix</b>					
			Internal hex	55020	10			