

SCALING

製品名	設定上限値
G1	0-7
G2	0-7
G3	0-7
G4	0-5
G5	0-6
G6	0-8
G8	0-7
G9	0-5
G11	0-6
G12	0-7
G13	0-6
G16	0-8

ENDODONTICS

製品名	設定上限値
E4	0-6
E4D	0-7
E7	0-6
E7D	0-6
E8	0-2
E8D	0-2
E15D	0-6
E11	0-4
E12	0-3

PERIO

製品名	設定上限値
P1	0-6
P1D	0-10
P2D	0-8
P3D	0-8
P10	0-10
P11R/L	0-10
P20	0-8
P21R/L	0-8
P25R/L	0-10
P26R/L	0-6
P40/P41	0-3
V-Tip V-P10 V-P11R/L V-P12 V-P26R/L	0-6

